GLUMA® Etch 35 Gel

(GB)	Instructions for use	2
(FR)	Mode d'emploi	6
(ES)	Instrucciones de uso	10
(GR)	Οδηγίες χρήσης	14
(CZ)	Návod k použití	18
(HU)	Használati utasítás	22
(LV)	Lietošanas instrukcija	26
(LT)	Naudojimo instrukcija	30
(PL)	Instrukcja obsługi	34
(BR)	Instruções de uso	38
(HR)	Upute za uporabu	42
(RO)	Instrucțiuni de utilizare	46
(TR)	Kullanma talimatı	50
(RU)	Инструкция по применению	54
(AE)	تعليمات الاستخدام	58
(CN)	使用说明	62
(KR)	사용 설명서	65
(SK)	Návod na použitie	68
(BG)	Инструкции за употреба	72
(TH)	คำแนะนำในการใช้	76
(UA)	Інструкція по застосуванню	80
(EE)	Kasutusjuhend	84
(SI)	Navodila za uporabo	88

KULZER
MITSUI CHEMICALS GROUP

Product description

GLUMA Etch Gel is an etching agent based on orthophosphoric acid. It is used to etch enamel and condition dentine in the context of adhesive filling therapy.

Composition

GLUMA Etch 35 Gel contains 35% phosphoric acid.

The product further contains thickening agents, pigments, water.

Indications

Enamel etching and enamel/dentine etching in the context of the total etch technique prior to the adhesive bonding of:

- Composite and compomer restorations
- Laboratory-fabricated, indirect restorations (e.g. inlays, crowns, bridges, veneers)
- Sealants

Contraindications

- Application to dentine close to the pulp (dentin layer less than 1 mm in thickness)
- Do not use in case of known hypersensitivity to GLUMA Etch

Side effects

This product or one of its components may in particular cases cause hypersensitive reactions. If suspected, information on the ingredients can be obtained from the manufacturer.

Interactions

Liner and base materials, such as calcium hydroxide, may be dissolved by the acid, depending on their type and composition.

Warnings/precautions

GLUMA Etch 35 Gel contains 35% phosphoric acid and can cause chemical burns. May irritate eyes. Avoid contact with gingiva, skin, eyes or neighbouring teeth. In case of contact with oral mucosa, rinse thoroughly with plenty of water. If swallowed, thoroughly rinse mouth with water and then drink plenty of water. Consult a physician. In case of contact with skin, thoroughly rinse with soap and water. In the event of contact with eyes, rinse thoroughly with water and seek medical advice. Wear suitable protective gloves. Avoid skin contact. Wear protective clothing. Wear safety glasses/face protection. Immediately remove clothing contaminated with GLUMA Etch 35 Gel.

Application

Prior to treatment, carefully clean tooth with fluoride-free polishing paste.

The use of a rubber dam is recommended.

1. After preparation in accordance with adhesive filling therapy, rinse the cavity and dry with oil-free air.
2. Starting from the bevelled enamel areas, apply GLUMA Etch to the entire cavity surface, including the dentine. Ensure that the enamel is etched for at least 20 sec and no more than 30 sec, and that the dentine is conditioned for at least 15 sec and no more than 20 sec. When using pay special

attention to complying with the dwell time on the dentine. Reducing the time to 15 sec in these areas is permissible to avoid overetching. If desired, only apply GLUMA Etch to the enamel and let set for a maximum of 30 sec. Etch unprepared enamel, e.g. in case of fissure sealing, for at least 30 sec and no more than 60 sec. In case of fluoridated enamel, etching for no more than 60 sec may also be required.

Caution: In case of deep cavities, we recommend covering the dentine near the pulp with a calcium hydroxide preparation and a base (e.g. glass ionomer cement) prior to etching.

3. Use a water sprayer to completely rinse off GLUMA Etch for 20 sec. The cavity must be completely free of any remnants.
4. Dry the cavity with oil-free air (do not dry out dentine!). The etched enamel surface should appear chalky white. In case of contamination of the tooth surface (e.g. with saliva), the etching must be repeated (for a maximum of 10 sec dwell time).
5. Perform further treatment with adhesive (e.g. GLUMA® 2Bond) and composite (e.g. Charisma® Composite) in accordance with manufacturer directions.

Instructions

GLUMA Etch Gel may be applied using a brush or cannula. Both brushes and cannulas are for one-time use only. To avoid cross infection, items intended for single use must be disposed of immediately after initial use. Release the plunger after applying the gel. Carefully reclose after removing product.

Special instructions

Only to be used by dentists and for its intended use.
Store out of reach of children.

Storage conditions

Store at temperatures below 25°C (77°F).

Do not use the material after its expiration date.

Please quote batch number and expiry date in all correspondence about the product.

Expiration date  and batch number : See note on syringe/packaging.

Only recycle completely empty packaging.

Remaining etching gel may be discharged with wastewater if highly diluted with water or neutralised.

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Kulzer GmbH