KAVITAN™ PRO High-Strength Glass-Ionomer Cement

en

Charakteristika:
KAVITAN PRO je a radiopaké chemically curing glass-ionomer (glass-polyalkenoat) cement, intended for semi-permanent fillings in posterior teeth. KAVITAN PRO features high strength, chemical adhesion to hard dental tissues and long-term fluoride ion release. Cement is highly viscous, non-sticky, easy to apply and to adapt (condense) in the cavity.

- Indications:**
• Semi-permanent fillings of Class I, II, III, not stressed by occlusal forces in contact points.
• Class V and III fillings in permanent teeth.
• Laminated fillings of Class I, II, III and V in permanent teeth - so-called sandwich technique.
• Fillings of Class I, II, III and V in deciduous teeth.
• Adhesive denture composite or amalgam fillings.
• Lost dentine build-ups prior to preparation for crowns, inlays or onlays or aesthetic veneers.
• Temporary fillings of Class I, II, III and V in permanent teeth.
• A.R.T. technique (Atraumatic Restorative Treatment).

Contraindications:
Do not use for Class IV fillings, carious lesions of Class II with conventional preparation cavity in permanent teeth and for cavities extending to occlusal tooth cusps. Do not use in cases of known hypersensitivity of the patient to any of the product components.

INSTRUCTION FOR USE:

Powder : liquid mixing ratio	3.5 : 1 (g:g)
Mixing time	30 - 45 sec
Working time (from the beginning of mixing to the initial of setting)	2 - 3 min
Manipulation time (from the end of mixing to the beginning of setting)	min. 1:30 min
Setting time (from the beginning of mixing to the initial setting)	min. 5 min

At 23 °C and relative humidity of 50 %.

Cement preparation:
• Mixing ratio: 3.5 : 1 (g : g), i.e. 1 measuring scoop of powder to 1 drop of liquid. Cement is liquid at 20 - 25 °C.

• Store the bottle with powder well before use. Do not compress powder in the measuring scoop. Overfill the scoop with the powder and level it off using the edge insert in the bottle neck.

• Hold the bottle vertically when dosing the liquor. You will get the optimum drop size when you squeeze the bottle gently. Any dried material must be removed from the dropper.

• Close the bottles immediately after closing of powder and liquid.

• Dose the powder and liquid on a disposable paper mixing pad or a glass slab (with a pre-dried stainless steel slab you can extend the working time).

• Preparation is incomplete if the powder is dry. Mix powder with liquid using a plastic or stainless steel spatula. Mixing time should be as short as possible (30-45 sec at most). Properly prepared cement is homogeneous, plastic and non-sticky.

Principles of clinical application:

• General preparation of the carious lesion, even creating macroretention (in the case of erosive or abrasive lesion no preparation is needed, only mechanical cleaning of the dentine surface is required).

• For repairing a pulpal wall thinner than 0.5 mm, use a Ca(OH)₂ based cement for lining.

• Adhesive will be enhanced by applying a conditioner - see Instructions for use for the product used.

• Apply the cement using a burnerish or a similar tool with high temperatures (or special applicators (Centrix tips, etc.). The material is non sticky, it can be shaped (condensed) using a ball pointed plugger of a suitable diameter).

• It is recommended to let the material set under a matrix.

• Applying, treat surface fillings with a protective lacquer (LC VARNISH).

• Remove excess material by manual tools (scapel, scalper, excavator) or rotating instruments.

• Final finishing, or etching for laminated fillings, if applicable, can start only after 15 minutes at the soonest.

Note:
Glass-ionomer cements are very susceptible to excessive moisture in the initial setting phase (disintegration in the mouth) and to overdrying during the same phase, therefore protective varnish should be always applied. The material tends to adhere to tools and it should be removed before it dries.

Warning:
Liquid: H319 Causes serious eye irritation. H315 Causes skin irritation. P262 Do not get in eyes, on skin, or on clothing. P270 Do not eat, drink or smoke when using this product.

Waste Disposal:
Dispose of unused powder and its packaging with ordinary worksite waste. Residues of unused liquid material should be transferred to a company authorized to dispose of dangerous waste materials.

Storage:
Store in a dry and dark place at 5 - 25 °C, in a well sealed container.

Shelf-life:
3 years

Indication on the package: Exp YYYY-MM (YYYY = year, MM = month).

Packaging:
15 g powder, shades A2, A3 + 10 g liquid, mixing pad, measuring scoop

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